

From **QR** **WOUNDSEAL™**

**STOP BLEEDING  
INSTANTLY!**

**Non-prescription Topical Powder**

***What is WoundSeal Powder?***

WoundSeal Powder is a non-prescription topical powder that controls external bleeding. It is composed of a hydrophilic polymer and potassium ferrate.

***How does it work?***

When the powder is poured onto a bleeding wound and manual pressure is applied, an “instant scab” or seal quickly forms over the wound.

***Powder + Pressure = Instant Scab***

**Instant Scab:**

- Stops Bleeding
- Seals and protects the wound
- Falls off naturally as the wound heals
- Is waterproof

***Safety Information:***

- WoundSeal Powder passed all safety and toxicity testing at an independent FDA licensed laboratory.
- Not OSHA-reportable as it is considered “first aid” treatment.
- Is not a drug; It is a non-prescription medical device.
- Is not biologically derived.



**WoundSeal Pour Pack**  
(was PRO QR for Cuts & Lacerations)  
For Easy-To-Reach Wounds -  
Cuts, Lacerations



**WoundSeal + Applicator**  
(was PRO QR for Nosebleeds)  
For Nosebleeds, Head Wounds and  
other Hard-To-Reach Wounds